

# Healthweigh® Indicators

## Load Cell Cable Connection

---

**Product:** Healthweigh Indicator

**Issue:** SAT appears on the indicator display

**Cause:** Improper cable connection

**Solution:** Connect the cable into the load cell connection port

---

When SAT appears on the display, the load cell cable may be plugged into the wrong connection port.




Manuals and additional resources are available from Rice Lake Weighing Systems at [www.ricelake.com/manuals](http://www.ricelake.com/manuals)

Warranty information is available at [www.ricelake.com/warranties](http://www.ricelake.com/warranties)

Use the following steps to ensure the cable is connected to the load cell connection port:


1. Remove the four screws securing the back cover to the indicator (see [Figure 1](#)).



*Figure 1. Back Cover Screw Locations*

2. Connect the load cell cable to the indicator by plugging it into the load cell connection port.
3. Secure the cable with the zip tie (see [Figure 2](#)).

 **IMPORTANT:** Do not plug the load cell cable into the RS-232 connection port on the right side.



*Figure 2. Indicator to Load Cell Cable Connection*

4. Ensure the cable routes through the notch in the top of the back cover (see [Figure 3](#)).



*Figure 3. Cable Through Back Cover Notch*

5. Reinstall the back cover of the indicator.



© Rice Lake Weighing Systems Content subject to change without notice.

230 W. Coleman St. • Rice Lake, WI 54868 • USA USA: 800-472-6703 • International: +1-715-234-9171