Medical Scales

TABLE OF CONTENTS

PLATFORM/FITNESS SCALE	212
MECHANICAL PHYSICIAN SCALES	213
MEDICAL INDICATOR, TITANIUM	214
150-10-7 Floor Level	215
150-10-5, 150-10-6	216
350-10-7, 350-10-8	217
550-10-1	218
350-10-4	219
350-10-2, 350-10-3	220
250-10-4	221
250-10-2	222
260-10-1	223
DIGITAL BABY SCALE	224
MECHANICAL CHAIR SCALE	225
PORTABLE BARIATRIC/WHEELCHAIR SCALE	226
PORTABLE BARIATRIC/WHEELCHAIR SCALE	228
BARIATRIC/WHEELCHAIR SCALE. Beam and Column Assembly	230
DIALYSIS SCALE	231
MEDICAL INDICATOR, PLATINUM	232



Platform/Fitness Scale





Legal for Trade, floor level Model 140-10-7N

Part#	Description	Price
125406	Load cell	Consult
131089	Indicator battery board (model 140-10-7N)	Consult
131090	Display/CPU board	Consult
102503	Indicator keypad	Consult
102506	Sealed lead acid battery	Consult
102507	AC adaptor, 110 VAC	Consult
131581	Load cell, NTEP Certified	Consult
106881	Indicator with tilt stand	Consult



Mechanical Physician Scales

Part#	Description	Price
109637	Height rod (RL-MPS, RL-MPS-10)	Consult
116380	Scale base foot (set of 5)	Consult
116381	Pillar screw with lock washer (set of 8)	Consult
175418	Black plastic cover for scale base	Consult
169569	Counterpoise set, one fine and one coarse weight	Consult




Model RL-MPS, RL-MPS-10

Medical Indicator, Titanium




Ref#	Part#	Description	Price
	131583	Indicator	Consult
10	164239	Front, 150 indicator kit	Consult
20	164240	Hardware	Consult
30	164241	Circuit	Consult
	204810	Circuit with Bluetooth 2.0	Consult
	204811	Circuit with BLE	Consult
40	164242	Back panel	Consult
50	164243	Table	Consult
60	164247	Back cover	Consult
70	164248	Cover	Consult
	133077	Power supply, A/C adapter	Consult



Ref#	Part#	Description	Price
10	164299	Feet (4)	Consult
20	164320	Platform without mat	Consult
30	164301	Mat	Consult
40	164321	Base	Consult
50	164304	Hardware, kit for load cell	Consult
60	164322	Cable clamp	Consult
71	164323	Load cell	Consult






150-10-5, 150-10-6




Ref#	Part#	Description	Price
10	164299	Feet	Consult
21	164300	Platform	Consult
30	164301	Mat	Consult
41	164303	Base	Consult
50	164304	Hardware	Consult
75	164306	Load cell	Consult
80	164305	Wheels	Consult
90	164307	Pillar	Consult
100	164308	Pillar, top	Consult
101	164309	Pillar, top	Consult
110	164310	Cable	Consult
130	164311	Display	Consult
140	164312	Display	Consult
150	164313	PCB	Consult
	207939	PCB with Bluetooth 2.0	Consult
	207548	PCB with BLE	Consult
160	164314	Back cover	Consult
170	164315	Battery	Consult
180	164316	Height rod	Consult
190	164317	Height rod	Consult
200	164318	Height rod	Consult
210	164319	Height rod	Consult



Ref#	Part#	Description	Price
10	164516	Feet	Consult
22	164518	Base	Consult
31	164521	Load cell	Consult
41	164534	Wheel	Consult
51	164578	Junction box	Consult
62	165029	Platform	Consult
66	165063	Platform	Consult
95	165072	Cable	Consult
110	165086	Spirit level assembly	Consult
130	165089	Carry handle	Consult
140	165090	Ramp	Consult








350-10-8

350-10-7




Ref#	Part#	Description	Price
2	164326	Castor	Consult
3	164327	Base	Consult
6	164328	Base	Consult
7	164329	Footrest	Consult
8	164330	Base	Consult
9	164331	Adapter	Consult
10	164332	Left arm	Consult
11	164333	Right arm	Consult
12	164334	Bracket	Consult
13	164324	Seat	Consult
14	164335	Handle	Consult
15	164336	Load cell and mounting screws	Consult
17	164337	Handles	Consult
18	164338	Cable	Consult



Ref#	Part#	Description	Price
10	164516	Feet	Consult
21	164518	Base	Consult
31	164521	Load cell	Consult
41	164534	Wheel	Consult
51	164578	Junction box	Consult
66.1	165064	Platform	Consult
72	165066	Joint	Consult
82	165069	Joint	Consult
	102636	Circular handle for hinge joint	Consult
92	165071	Cable	Consult
104	165085	Handrail kit	Consult
110	165086	Spirit level assembly	Consult
121	165088	Seat kit	Consult
130	165089	Carry handle	Consult
140	165090	Ramp	Consult
150	165091	Rubber stop	Consult
160	165092	Back seat	Consult






350-10-2, 350-10-3





Ref#	Part#	Description	Price
10	164516	Feet	Consult
21	164518	Base	Consult
31	164521	Load cell	Consult
41	164534	Wheel	Consult
51	164578	Junction box	Consult
62.1	165030	Platform	Consult
66.1	165064	Platform	Consult
71	165066	Right joint	Consult
81	165069	Left joint	Consult
	102636	Circular handle for hinge joint	Consult
92	165071	Cable	Consult
102	165084	Handrail kit	Consult
110	165086	Spirit level assembly	Consult
130	165089	Carry handle	Consult
140	165090	Ramp	Consult
150	165091	Rubber stop	Consult




Ref#	Part#	Description	Price
10	164516	Feet	Consult
20	164517	Base	Consult
31	164521	Load cell	Consult
40	164527	Wheel	Consult
51	164578	Junction box	Consult
60	165028	Platform	Consult
70	165065	Joint	Consult
80	165068	Joint	Consult
92	165071	Cable	Consult
101	165081	Handrail kit	Consult
110	165086	Spirit level assembly	Consult
120	165087	Seat kit	Consult









Ref#	Part#	Description	Price
	195888	Wheels (set of 2)	Consult



Digital Baby Scale



Model RL-DBS/RL-DBS-2

Part#	Description	Price
112692	120 VAC adaptor	Consult
112695	230 VAC adaptor	Consult
121333	Battery cover	Consult


Mechanical Chair Scale

Part#	Description	Price
116392	Front wheel, caster (set of 2)	Consult
116393	Back wheel, caster (set of 2)	Consult
116394	Left arm	Consult
116395	Right arm	Consult



Portable Bariatric/Wheelchair Scale

Dual Ramp



Portable Bariatric/Wheelchair Scale Dual Ramp Model RL-350-6


Part Number/Price


Ref#	Part#	Description	Price
	100719	RS-232 serial cable	Consult
3	105373	Wheel kit (includes wheel, torque nut, screw, scotch mate hoops)	Consult
4	106953	Cable assembly	Consult
5	113030	Ramp	Consult
6	113031	Junction box cover	Consult
7	113033	Handle insert	Consult
8	114614	Foot	Consult
9	14626	Nut, kep 8-32NC hex	Consult
10	14839	Screw, mach 6-32NC x 1/4	Consult
11	14857	Screw, mach 8-32NC x 1/4	Consult
12	14862	Screw, mach 8-32NC x 3/8	Consult
14	58983	Cord grip	Consult
15	84838	Load cell	Consult
16	88958	Junction box board	Consult

Note: See indicator parts on page 214

(See Diagram Portable Bariatric/Wheelchair Scale Dual Ramp)

Portable Bariatric/Wheelchair Scale Dual Ramp






Portable Bariatric/Wheelchair Scale Dual Ramp

MEDICAL SCALES

Portable Bariatric/Wheelchair Scale

Single Ramp



Portable Bariatric/Wheelchair Scale Single Ramp with Wheel Stop Model RL-350-5


Part Number/Price

Ref#	Part#	Description	Price
	100719	RS-232 serial cable	Consult
3	105373	Wheel kit (includes wheel, torque nut, screw, scotch mate hoops)	Consult
4	106953	Cable assembly	Consult
5	113030	Ramp	Consult
6	113031	Junction box cover	Consult
7	113033	Handle insert	Consult
8	114413	Wheelchair stop	Consult
9	114614	Foot	Consult
10	14626	Nut, kep 8-32NC hex	Consult
11	14635	Nut, lock 1/4-20NC hex	Consult
12	14839	Screw, mach 6-32NC x 1/4	Consult
13	14857	Screw, mach 8-32NC x 1/4	Consult
14	14862	Screw, mach 8-32NC x 3/8	Consult
15	14955	Screw, cap 1/4-20NC x 1/2	Consult
17	58983	Cord grip	Consult
18	84838	Load cell	Consult
19	88958	Junction box board	Consult

Note: See page 214 for indicator parts

(See Diagram Portable Bariatric/Wheelchair Scale Single Ramp)

Portable Bariatric/Wheelchair Scale Single Ramp



Portable Bariatric/Wheelchair Scale Single Ramp


Bariatric/Wheelchair Scale Beam and Column Assembly



Ref#	Part#	Description	Price
1	18756	Foot	Consult
2	21507	Load cell	Consult
3	30532	Load cell bolts	Consult
4	73428	Junction box	Consult
	106953	Cable, junction box to indicator	Consult
	100719	RS-232 serial cable	Consult

Note: See page 214 for indicator parts (See Diagram Wheelchair-F)

Bariatric/Wheelchair scale Model SB-1150



Wheelchair-F


Dialysis Scale

Ref#	Part#	Description	Price
1	18756	Foot	Consult
2	21507	Load cell	Consult
3	30532	Load cell bolts	Consult
4	73428	Junction box	Consult
5	107901	Handrail tube and base plate	Consult
6	107902	Washer	Consult
7	15019	Hex head screw	Consult
	106953	Cable, junction box to indicator	Consult
	100719	RS-232 serial cable	Consult





Dialysis Scale Model SD-1150-WP



Dialysis-A

Medical Indicator, Platinum



Part Number/Price

Ref#	Part#	Description	Price
2	116382	Tilt stand	Consult
3	116384	Tilt stand washer	Consult
6	116383	Tilt stand knob	Consult
9	102519	Pan head screw	Consult
13	102519	Pan head screw	Consult
14	105378	Indicator house	Consult
17	102506	Sealed lead acid battery	Consult
19	102505	Indicator battery board	Consult
19	102505	Indicator battery board (model 140-10-7N)	Consult
21	102504	Display/CPU board	Consult
25	116385	Indicator lens	Consult
26	105379	Indicator front cover	Consult
27	102503	Indicator keypad	Consult
	114582	Jumper	Consult

(See Diagram Indicator-A)

