

1351-10-2

IEELCHAIR SCALE WITH HANDRAIL

The Rice Lake Healthweigh® H351-10-2 single-ramp wheelchair scale with handrail allows high capacity, accurate weighing for a wide range of patients. Convenient keypad tare reduces errors and increases efficiency. Foldable and lightweight, this sturdy scale is easily transported on two heavy-duty rear wheels to any location.

FEATURES

- Two wheels for portability
- Single ramp
- Nonslip platform
- Large LCD display
- Motiontrap[™] movement compensation technology
- Unit of measure kg or lb

- Serial and USB data communication ports
- Tare function
- Hold function
- BMI (body mass index) function
- Low battery annunciator
- Easy to store, folded depth is 30 cm



Platform Dimensions Without Ramps:

 $81.3 \times 74.9 \times 6.1 \text{ cm}$

Rechargeable NiMH battery pack with 9 VDC adaptor, 120V/230V auto switching

CEE 7/16 Plug (0 0)

Battery Use:

25 hours of continuous use with batteries Low battery annunciator

LED indicates batteries are charging (red) or fully charged (green)

Automatic power-off can be configured

Operating Temperature:

10 °C to 40 °C

Display:

5-digit LCD display, 1.9 cm digit height

Approvals:

CE

c-Met, E113986

OIML Class III

Warranty:

Two-year limited





WHEELCHAIR SCALE


MODEL	CAPACITY	PART#
H351-10-2	$300 \times 0.1 \text{kg}$	217646

SHIPPING WEIGHT (WITHOUT PALLET)

MODEL#	PART#	L (cm)	W (cm)	H (cm)	PACKAGE WGT. (kg)	DIM WGT. (kg)	SHIPPING WGT. (kg)
H351-10-2	217646	107	104	36	41	79	79

DIMENSIONS - cm

A	74.9	D	23.5
В	81.3	E	95.3
C	6.1	F	114.3



An ISO 9001 Registered Company ©2024 Rice Lake Weighing Systems 01/24 Specifications subject to change without notice. Visit ricelake.com