H351-10-3 Dual-ramp wheelchair scale with handrail

The Rice Lake Healthweigh[®] H351-10-3 dual-ramp wheelchair scale with handrail allows high capacity, accurate weighing for a wide range of patients. Convenient keypad tare reduces errors and increases efficiency. Foldable and lightweight, this sturdy scale is easily transported on two heavy-duty rear wheels to any location.



- Two wheels for portability
- Dual ramp
- Nonslip platform
- Large LCD display
- Motiontrap[™] movement compensation technology
- Unit of measure kg or lb

SPECIFICATIONS

Capacity: 300 x 0.1 kg

Platform Dimensions Without Ramps: (W × L × H) $81.3 \times 74.9 \times 6.1$ cm

Power:

Rechargeable NiMH battery pack with 9 VDC adaptor, 120V/230V auto switching CEE 7/16 Plug

- Serial and USB data communication ports
- Tare function
- Hold function
- BMI (body mass index) function
- Low battery annunciator
- Easy to store, folded depth is 30 cm

Battery Use:

25 hours of continuous use with batteries Low battery annunciator LED indicates batteries are charging (red) or fully charged (green) Automatic power-off can be configured

Operating Temperature: 10 °C to 40 °C

Display: 5-digit LCD display, 1.9 cm digit height

Approvals: CE c-Met, E113986 OIML Class III

Warranty: Two-year limited



WHEELCHAIR SCALES



WHEELCHAIR SCALE

MODEL	CAPACITY	PART#
H351-10-3DI	300 × 0.1 kg	217571

SHIPPING WEIGHT

MODEL #	PART #	L (cm)	W (cm)	H (cm)	PACKAGE WGT. (kg)	DIM WGT. (kg)	SHIPPING WGT. (kg)
H351-10-3	217571	130	107	36	46	98	98

DIMENSIONS - cm

Α	123.2	D	23.6
В	81.3	E	96.5
С	6.1	F	114.3



An ISO 9001 Registered Company ©2024 Rice Lake Weighing Systems 01/24 Specifications subject to change without notice. Visit ricelake.com