H351-10-8 Dual-ramp wheelchair platform scale

The Rice Lake Healthweigh® H351-10-8 dual-ramp wheelchair platform scale allows high capacity, accurate weighing for a wide range of patients. Convenient keypad tare reduces errors and increases efficiency. This sturdy scale is easily transported on two heavy-duty rear wheels to any location.

FEATURES

- Two wheels for portability
- Dual ramp
- Nonslip platform
- Large LCD display
- Motiontrap[™] movement compensation technology
- Unit of measure kg or lb

SPECIFICATIONS

Capacity: 300 x 0.1 kg

Platform Dimensions Without Ramps: (L × W × H) 74.9 × 81.3 × 6.1 cm

Power:

Rechargeable NiMH battery pack with 9 VDC adaptor, 120V/230V auto switching CEE 7/16 Plug



- Serial and USB data communication ports
- Tare function
- Hold function
- BMI (body mass index) function
- Low battery annunciator

Battery Use:

25 hours of continuous use with batteries Low battery annunciator LED indicates batteries are charging (red) or fully charged (green) Automatic power-off can be configured

Cable Length: 2.1 m cable between scale base and indicator

Operating Temperature: 10 °C to 40 °C

Display: 5-digit LCD display, 1.9 cm digit height

Approvals:

CE c-Met, E113986 OIML Class III

Warranty: Two-year limited

WHEELCHAIR SCALES



WHEELCHAIR SCALE

MODEL	CAPACITY		PART#				
H351-10-8	300 × 0.1 kg		217663				
SHIPPING WEIGHT							
		DAOKAOE					

MODEL #	PART #	L (cm)	W (cm)	H (cm)	WGT. (kg)	(kg)	WGT. (kg)
H351-10-8		100		10	00	55	55

DIMENSIONS - cm

Α	123.2	D	23.6
В	81.3	E	96.5
C	6.1		

Rear Wheels D

An ISO 9001 Registered Company ©2024 Rice Lake Weighing Systems 01/24 Specifications subject to change without notice. Visit ricelake.com