

# MASS REFERENCE STANDARDS / METRIC SETS

## ESSENTIAL WEIGHT SETS

### Class 0

Set Range	Part #	Price
2 g, 200 mg, 1 mg	<b>186177</b>	Consult
5 g, 500 mg, 1 mg	<b>186178</b>	Consult

### Class 1

Set Range	Part #	Price
20 g, 2 g, 20 mg	<b>186179</b>	Consult
100 g, 10 g, 20 mg	<b>186180</b>	Consult
100 g, 10 g, 200 mg	<b>186181</b>	Consult
200 g, 20 g, 20 mg	<b>186182</b>	Consult
200 g, 20 g, 200 mg	<b>186183</b>	Consult
200 g, 20 g, 2 g	<b>186184</b>	Consult
500 g, 50 g, 200 mg	<b>186185</b>	Consult
500 g, 50 g, 2 g	<b>186186</b>	Consult
1 kg, 100 g, 2 g	<b>186187</b>	Consult
1 kg, 100 g, 20 g	<b>186188</b>	Consult
2 kg, 200 g, 2 g	<b>186189</b>	Consult
2 kg, 200 g, 20 g	<b>186190</b>	Consult
5 kg, 500 g, 2 g	<b>186191</b>	Consult
5 kg, 500 g, 20 g	<b>186192</b>	Consult
100 g, 10 g, 2 g	<b>191089</b>	Consult
500 g, 50 g, 20 mg	<b>191091</b>	Consult
50 g, 5 g, 200 mg	<b>191092</b>	Consult
50 g, 5 g, 2 g	<b>191093</b>	Consult

### Standard Features

- Polished stainless steel Type II design
- Sealing cavity
- Density is approximately 7.95 g/cm<sup>3</sup>
- Special PVC weight case
- Precision tweezers included
- Configurations are ideal for calibration, shift tests and sensitivity tests
- Other configurations available by special order
- Meets USP-41 requirements for intermediate checks
- Includes ISO 17025 Accredited Certificate of Calibration



200g PVC Set Shown, Other Sets Will Vary

### Accessories

Description	Part #	Price
Weight glove, white cotton	<b>11714</b>	Consult
Tweezer	<b>11711</b>	Consult
Premium tweezer	<b>152414</b>	Consult
Lifter bar	<b>11747</b>	Consult

ASTM  
PRECISION

OIML  
PRECISION

BALANCE  
CALIBRATION

STAINLESS  
FIELD WEIGHTS

CAST IRON  
FIELD WEIGHTS

SPECIALTY

ACCESSORIES