

NATIONAL TYPE EVALUATION PROGRAM

Certificate of Conformance for Weighing and Measuring Devices

For:

Indicating Element Digital Electronic

Model: SCT20-XX (See below)

n_{max}: 5000

Accuracy Class: III

Submitted By:

Rice Lake Weighing Systems

230 W Coleman St.

Rice Lake, WI 54868

Tel: 715-234-9171 x 5322

Fax: 715-234-6967

Contact: Paul A. Lewis, Sr. Email: plewis@ricelake.com

Standard Features and Options

Automatic Zero Tracking (AZT) Semi-Automatic (Push Button) Zero Programmable Tare Semi-Automatic (Push Button) Tare DC Power Supply 12-24 VDC Gross/Net Display Linearity Calibration Points (5)

Remote Calibration

LED Display

XX Suffix Designates Communication Options:

SCT20-AN (analog + serial RS485)

SCT20-485 (serial RS485)

SCT20-PB (PROFIBUS)

SCT20-PN (PROFINET/IO)

SCT20-CO (CANOPEN)

SCT20-DN (DEVICENET)

SCT20-IP (ETHERNET/IP)

SCT20-EC (ETHERCAT)

SCT20-TCP/IP (ETHERNET/TCP/IP)

SCT20-MB (MODBUS/TCP)

SCT20-CC (CC-LINK)

SCT20-SERCOS (SERCOS-III)

SCT20-PL (POWERLINK)

Temperature Range: -10 °C to 40 °C (14 °F to 104 °F)

This device was evaluated under the National Type Evaluation Program and was found to comply with the applicable technical requirements of "NIST Handbook 44: Specifications, Tolerances and Other Technical Requirements for Weighing and Measuring Devices." Evaluation results and device characteristics necessary for inspection and use in commerce are on the following pages.

Ronald Haves

Chairman, NCWM, Inc.

John Gaccione Committee Chair, National Type Hvaluation Program Committee

Issued: September 23, 2014

1135 M Street, Suite 110 / Lincoln, Nebraska 68508

The National Conference on Weights and Measures (NCWM) does not approve, recommend or endorse any proprietary product or material, either as a single item or as a class or group. Results shall not be used in advertising or sales promotion to indicate explicit or implicit endorsement of the product or material by the NCWM.





Rice Lake Weighing Systems

Indicating Element / SCT20-XX

Application: A general purpose indicating element to be interfaced with an NTEP certified and compatible weighing element.

<u>Identification</u>: The adhesive identification badge is located on the top of the indicator. If the badge is removed the word "VOID" repeatedly appears.

<u>Sealing</u>: The indicator is sealed by using a self-destructive adhesive label place over the seam between the top and bottom half of the indicator. The seal prevents the top and bottom half from being separated allowing access to the jumper shown below from being moved. When the jumper is on, the indicator can be calibrated and configured. The jumper is off the indicator cannot be calibrated or configured. Also when in calibration or configuration mode, the annunciators continually scroll to show the indicator is not in the regular weighing mode. See pictures below for sealing method and jumper location.



Jumper Location


<u>Test Conditions</u>: This certificate was issued based upon the following tests and upon information provided by the manufacturer. The emphasis of this evaluation was on the device design, operation, marking requirements, performance, and compliance with influence factors. The indicator was interfaced with a Rice Lake BM 1212 (NTEP CC 95-072) weighing/load receiving element to verify zero, discrimination zone of uncertainty, and motion detection requirements. A load cell simulator was used to perform several increasing/decreasing tests, DC voltage testing at 12 VDC and 24 VDC, and temperature testing over a range of -10 °C to 40 °C (14 °F to 104 °F).

Evaluated By: T. Buck (OH)

Type Evaluation Criteria Used: NIST Handbook 44 Specifications, Tolerances, and Other Technical Requirements for Weighing and Measuring Devices, 2014 Edition. NCWM Publication 14 Weighing Devices, 2014 Edition.

Conclusion: The results of the evaluation and information provided by the manufacturer indicate the device complies with applicable requirements.

Information Reviewed By: J. Truex (NCWM)







Rice Lake Weighing Systems Indicating Element / SCT20-XX

Example of Device:

