

Mounting Assembly Component Identification







Conventional Tension Mounting Assembly

Our popular Isolated Tension Cell Mount (ITCM) assembly and other conventional suspended weighing system assemblies employ a variety of serviceable hardware components. Rice Lake Weighing Systems offers a comprehensive selection of these items for all of your suspended vessel service needs. These drawings are designed to enable you to quickly identify and differentiate components from both of these common designs. All corresponding replacement hardware is located on the following pages.

Jam Nuts

Thread Size	Part#	Price
Mild Steel, Fine Thread		
1/4-28 UNF	14639	Consult
3/8-24 UNF	14651	Consult
1/2-20 UNF	14664	Consult
3/4-16 UNF	14687	Consult
1-12 UNF	14702	Consult
1-14 UNF	14703	Consult
1 1/4-12 UNF	14708	Consult
1 1/2-12 UNF	14711	Consult
Mild Steel, Coarse Thread		
1/4-20 UNC	14636	Consult
3/8-16 UNC	14649	Consult
1/2-13 UNC	14661	Consult
5/8-11 UNC	14676	Consult
3/4-10 UNC	14686	Consult
1-8 UNC	14701	Consult
1-1/4-7 UNC	14707	Consult

S	Standard Features	
•	Coarse or fine thread	
_	7inc plated mild stool or stai	

- Zinc-plated mild steel or stainless steel
- **Thread Size** Part# Price Stainless Steel, Fine Thread 1/4-28 UNF 14640 Consult 1/2-20 UNF 14665 Consult 3/4-16 UNF 14688 Consult Stainless Steel, Coarse Thread 3/8-16 UNC 14653 Consult 1/2-13 UNC 14662 Consult 5/8-11 UNC 14680 Consult 3/4-10 UNC 14690 Consult